wellion® WAVE

HANDBUCH

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OWNERS GUIDE

TABLE OF CONTENTS

INTRODUCTION	23
ABOUT BLOOD PRESSURE	24
What is blood pressure?	24
What is hypertension and how is it controlled?	24
Why measure blood pressure at home?	24
WHO blood pressure classification	24
Blood pressure variations	24
NOTES ON SAFETY	25
Explanation of the icons	25
PRECAUTIONS BEFORE USE	26
PARTS IDENTIFICATION	27
INSERT OR REPLACE BATTERIES	28
TIME SETUP	29
UNIT CONVERSION mmHg/kPa DISPLAY	29
WHO BLOOD PRESSURE CLASSIFICATION DISPLAY	30
ATTACHING THE WRIST CUFF	30
Fastening the wrist cuff	30
How to take proper measurements	30
HOW TO MEASURE BLOOD PRESSURE	31
Read Memory	31
Delete Memory	31
CLEANING AND MAINTENANCE	32

TROUBLESHOOTING	33
SPECIFICATIONS	34
STATEMENT	34
FEATURE OF THE PRODUCT	35
REQUESTS FROM MANUFACTURER	35
REGARDING EMC STATEMENT	36
MANUFACTURER'S WARRANTY	36
Appendix 1 Guidance and Manufacturer Declaration Tables	37
CALIBRATION METHOD	41

INTRODUCTION

The Wellion WAVE Blood Pressure Monitor uses the oscillometric method of blood pressure measurement.

The electronic Wellion WAVE Blood Pressure Monitor for automatic measurement is intended for use by medical professionals or at home to monitor and display diastolic and systolic blood pressure and pulse rate, with an air wrist cuff buckled around the wrist according to the instructions in the chapter "ATTACHING THE WRIST CUFF."

The expected lifespan of the product is 5 years.

The product complies with the electromagnetic compatibility requirement of IEC 60601-1-2, safety standards of IEC 60601-1 and performance requirements of IEC 80601-2-30 as specified in Regulation (EU)2017/745.

ABOUT BLOOD PRESSURE

What is blood pressure?

Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands.

Blood pressure is measured in millimetres of mercury (mmHg) or in kilopascals (kPa).

What is hypertension and how is it controlled?

Hypertension is an abnormally high arterial blood pressure. If this is left unattended it can cause many health problems including stroke and heart attack. Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision.

To prevent hypertension or to keep it under control:

- Do not smoke
- Reduce salt and fat intake
- Maintain proper weight
- Exercise regularly
- Have regular physical check-ups

Why measure blood pressure at home?

Blood pressure measured at a clinic or doctor's office may cause apprehension and produce an elevated reading, 25 to 30 mmHg higher than the value measured at home. Home measurement reduces the effects of outside influences on blood pressure readings, supplements the doctor's readings and provides a more accurate, complete blood pressure history.

WHO blood pressure classification

Standards for assessment of blood pressure, without regard to age, have been established by the World Health Organization (WHO) as shown (Graph taken from Journal of Hypertension 1999, Vol 17 No.2). It is recommended to take blood pressure measurements at approximately the same time each day, e.g. in the morning before breakfast.



Blood pressure variations

It is common for blood pressure readings to vary throughout the day. Normal fluctuations in blood pressure are typically caused by factors such as physical activity, stress levels, and even changes in weather.

Generally, blood pressure is highest during the day and lowest at night.

Take measurements at the same time every day using the procedure described in this manual, and know your normal blood pressure. Many readings give a more comprehensive blood pressure history.

Be sure to note date and time when recording your blood pressure. Consult your doctor to interpret your blood pressure data.

NOTES ON SAFETY

Please find below the meaning of some key signs and symbols providing necessary information to ensure the end user's safety as well as longevity of the device.

EN

Explanation of the icons

	The icon indicates prohibitions (what you should not do). Specific actions prohibited are indicated either by text or a picture within or near the circle. The left icon refers to "general prohibition".
	The icon indicates something that is compulsory (what must always be observed). Compulsory actions are indicated either by text or pictures within or near the circle. The left icon refers to "general compulsion".
	Do not disassemble
	Type BF applied part Electrically connected to patient but not directly to heart
	Caution
	Contact your local authorities to determine the proper method for disposal of potentially bio hazardous parts and accessories.

	Indicates a medical device that needs to be protected from moisture.
	Please refer to the instructions for use
	Model number
IP 20	Ingress protection degree as defined by IEC 60529 IP20: Product is touchproof and will be resistant to dust or objects that are over 12mm in size (e.g. persons fingers), but is not waterproof.
	MR Unsafe An item marked MR unsafe is known to pose hazards in all MR environments.
	Consult instructions for use
	Medical device


Only a healthcare professional is qualified to interpret blood pressure measurements. Self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. The device should not be used to decide on first aid measures nor for continuous monitoring of blood pressure. This device is not suitable for measurement in any case of medical emergency or during operations. Although portable, it should not be operated in transit (e.g. patient transport). It is intended for use at home or in any stationary location. Please press the "ON/OFF" button to stop the device from working when you feel uncomfortable with the wrist, or if the air is inflating abnormally without stop.	 Caution 
The device should not be operated by children under 18 or people with restricted physical, sensory, or mental acumen so to avoid hazardous situations.	
Do not use the device for any purpose other than measuring blood pressure:	
During measurement, make sure that no mobile phone or any other electrical devices emitting electromagnetic fields is in the surrounding area. This may result in incorrect operation of the device and/or cause an inaccurate reading.	
The device is prohibited from being used during movement.	
Do not use the equipment outdoors or in wet areas.	
Do not disassemble, repair, or remodel the main unit or the cuff of the blood pressure monitor. This will lead to malfunctioning of the device.	

PRECAUTIONS BEFORE USE


1. If you are taking medication, consult your doctor to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting your doctor.
2. For people with irregular or unstable peripheral circulation problems due to diabetes, liver disease, hardening of the arteries, etc., there may be differences in blood pressure values measured at the upper arm versus at the wrist.
3. Measurements may be impaired if this device is used near televisions, microwave ovens, X-ray, mobile phone equipment or other devices with strong electrical fields. To prevent such interference, use the monitor at a sufficient distance from such devices or turn them off.
4. Wash your hands before use.

5. Do not measure on the same arm which uses simultaneously another medical monitoring ME Equipment; otherwise it could cause loss of function.
6. Consult your doctor if unexpected readings are obtained, also please refer to "Trouble shooting" of the manual.
7. As people generally are more relaxed at home than at the doctor's office, blood pressure may be lower at home.
8. Cuff pressure range: 0-299mmHg


PARTS IDENTIFICATION

Accessory


Case


Manual

Symbols on display

INSERT OR REPLACE BATTERIES

1. Remove the battery cover.
2. Insert new batteries into the battery compartment as shown, making sure to match the correct polarity.
3. Close the battery cover. Use only LR03, AAA batteries.


Empty batteries must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Caution:

- Insert the batteries as shown in the battery compartment. If not, the device will not work.
- When „LOW BATTERY“ symbol  blinks in the display, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery lifespan, or cause the device to malfunction.
„LOW BATTERY“ symbol  does not appear when the batteries run out.
- Please ensure to match the correct polarity when replacing batteries.
- Batteries, which have a fluid on the surface or have been modified, must not be inserted into the product.
- Battery short circuit must be prevented.
- Battery lifespan varies with the ambient temperature and may be shorter at low temperatures.
- The batteries may leak and cause a malfunction.
- Use only batteries corresponding to given specifications. The batteries provided with the device are for testing monitor performance and may have a shorter lifespan.
- Used batteries may leak and damage the main unit. Please observe the following points:
 - Remove batteries if the device is not likely to be used for some time (approximately three months or more).

TIME SETUP

1. Press the „SET“ button to enter the time setting menu.
2. To enter the time setting mode when switched off, press and hold the “SET” button until the year flashes on the LCD.
3. Press the “MEM” button to adjust the year, then press the “SET” button again to save your settings and enter the month setting mode.
4. Press the “MEM” button to adjust the month. Repeat the same steps to adjust day/hour/minute until setting is completed.



UNIT CONVERSION mmHg/kPa DISPLAY

Blood pressure readings can be displayed in either one of the following units: mmHg or kPa (mmHg being the preset unit)

To enter the unit settings menu press the „ON/OFF“ button for 10 sec. Choose the unit of your liking by pressing the „MEM“ button and then exit the menu by pressing the „ON/OFF“ button again. The current unit is displayed on the LCD screen.

Please note that all memorized results change their unit as well.

WHO BLOOD PRESSURE CLASSIFICATION DISPLAY


Diastolic blood pressure Reference material:
journal of hypertension 1999, vol 17 No.2

■	—	Grade 3 hypertension (severe)
■	—	Grade 2 hypertension (moderate)
■	—	Grade 1 hypertension (mild)
■	—	High-normal
■	—	Normal
■	—	Optimal

ATTACHING THE WRIST CUFF

Fastening the wrist cuff

1. Wrap the wrist cuff around your wrist about 1-2 cm below your hand as shown in the figure below.
2. Fasten the wrist cuff tightly by using the Velcro strip. For proper measurements, fasten the wrist cuff tightly and measure on a bare wrist.

How to take proper measurements


For best accuracy in blood pressure measurement:

- Sit comfortably at a table. Rest your wrist on the table.
- Relax for about 5 to 10 minutes before making a measurement.
- Raise your hand so that the wrist cuff is at the same level as your heart.
- Remain still and keep quiet during measurement.
- Do not measure right after physical exercise or a bath.
- Measure your blood pressure at about the same time every day.

HOW TO MEASURE BLOOD PRESSURE

1. Fasten the wrist cuff according to the instructions in "ATTACHING THE WRIST CUFF".
2. Press the "ON/OFF" button; all icons automatically appear for two seconds on the display, then the device switches to measurement and either displays „0“ or shows the last measurement record.
3. When the measurement starts, the cuff in the strap will automatically inflate. The heartbeat icon will flash on the display. When the measurement is completed, the results will be displayed.


EN

Read Memory

Press the "MEM" button while the meter is off to retrieve saved data. The meter shows the average value of the latest 3 measurements on the display, indicated by the „**AUS**“ symbol in the right upper corner of the display. In order to retrieve your individual blood pressure values, press the „MEM“ button again. To go through the memory, press „MEM“ for navigating up or „SET“ for navigating down. It is also possible to review your measurements directly after measurement by pressing the „MEM“ button.


Delete Memory

Enter the "READ MEMORY" mode, press and hold the "MEM" button for five seconds until the LCD displays "no". This indicates that all stored results have been deleted successfully.

CLEANING AND MAINTENANCE

- Keep the monitor in the storage case when not in use.
- Do not fold the wrist cuff too tightly.
- Clean the monitor with a dry, soft cloth. Do not use any cleaning solution.
- Do not submerge the device or any components in water.
- Store the device and the components in a clean and safe location.
- The cleaning steps for the cuff are as follows:
 - Completely wipe the inner side of the cuff (the side that contacts skins) with a soft cloth lightly moistened with 75% Ethyl alcohol 3 times. Replace the soft cloth after each wipe.
 - Then air dry the cuff.

Caution:

- Do not subject the monitor to extreme hot or cold temperatures, humidity or direct sunlight.
- Do not drop this monitor or subject it to strong impact.
- Remove the batteries if the unit will not be used for three months or longer.
- Always replace the whole lot of batteries at the same time.

Information:

- The device was designed and manufactured for a long service life. However, it is generally recommended to have the device inspected and calibrated every 2 years, to ensure proper function and performance.

TROUBLESHOOTING

If you have trouble in using the product please check the following points first.

EN

Error	Possible cause	How to correct
Device doesn't power on when ON/OFF button is pressed or the battery icon  flashes	No batteries installed.	Insert batteries
	Batteries completely drained.	Replace batteries
	Polarity of batteries is not matched correctly.	Make sure to match the polarity of the batteries correctly.
E1: cuff cannot be inflated properly	Check your wrist cuff for any air leakages.	Replace your wrist cuff with a new one
E2: Movement during measurement.	Motion, trembling and shivering affects the measurement reading.	Remain still and do not talk during measurement.
E3: inflating pressure is too high	Pressure value of more than 299 mmHg.	Start a re-measurement or return the device for recalibration.
E4: unusual fluctuation of blood pressure values	Motion, trembling and shivering affects the measurement reading.	Remain still and do not talk during measurement.
Battery icon  displayed	Batteries have low power.	Replace batteries and repeat the measurement
The systolic or diastolic pressure value is too high	The wrist cuff was held lower than your heart.	Keep correct position and gesture to measure again
	The wrist cuff was not attached properly.	
	You moved your body or spoke during measurement.	
The systolic or diastolic pressure value is too low	The wrist cuff was held higher than your heart.	
	You moved your body or spoke during measurement.	

SPECIFICATIONS

Measuring Method:	Oscillometric Measurement	
Measuring Range:	Pressure: 30 - 280 mmHg	Pulse: 40 - 199 beats/minute
Memory:	Stores up to 90 blood pressure readings	
Operating condition:	+ 5°C - + 40°C; 15% RH - 93% RH Atmospheric pressure: 70kPa - 106kPa	
Storage condition:	- 20°C - + 55°C; 0% RH - 93% RH Atmospheric pressure: 50kPa - 106kPa	
Classification:	Type BF	
Indication:	Digital LCD display	
Accuracy:	Static Pressure: ± 3 mmHg	Pulse: $\pm 5\%$
Power supply:	2x 1,5V Batteries (LR03 or AAA), approximately 200 measurements (using new alkaline batteries)	
Dimensions:	Approx.: 72(W) X 67(H) X 33(D) mm	
Weight:	Approx.: 130g, excluding batteries	
Wrist circumference:	13,5 - 19,5 cm	

1. Type of protection against electric shock: INTERNALLY POWERED EQUIPMENT.
2. Degree or protection against electric shock: TYPE BF APPLIED PART.
3. Mode of operation: CONTINUOUS OPERATION.
4. Equipment not suitable for category AP&APG equipment use in presence.

STATEMENT

The system might not meet its performance specifications if stored or used outside the temperature and humidity as mentioned below:

- Operating conditions: +5°C ~ +40°C; 15%RH ~ 93%RH; 70kPa ~ 106kPa
- Storage conditions: -20°C ~ +55°C; 0%RH ~ 93%RH

FEATURE OF THE PRODUCT

1. Stores up to 90 blood pressure readings.
2. Large and clear LCD display.
3. WHO blood pressure classification display.
4. Easy to use: Press a single button to start automatic measurement; Both measurement values and measurement time are recorded.
5. Automatically turns off (within 1 minute) to save power.

REQUESTS FROM MANUFACTURER

- Do not take blood pressure measurements more often than three times in a row. In between measurements there should be a resting interval of at least 5 minutes in order to prevent blood flow interference.
- DO NOT take blood pressure measurements more often than necessary, do not exceed a number of 6 measurements per day.
- Do not apply the cuff over a wound as this can cause further injury.
- Do not measure on the wrist which is on the side of a mastectomy, otherwise it could cause injury.
- Observe the air pressure value from the LCD display.
- When measuring, the device cannot exceed 280 mmHg, otherwise please press the "ON/OFF" button to stop.
- Do not use force to bend the wrist cuff.
- DO NOT drop or subject the device to severe shocks or vibrations.
- Please use only accessories and detachable parts specified/ authorised by the manufacturer. Otherwise, it may cause damage to the unit or pose danger to the user/patients.
- For service information, parts list etc., please contact the dealer.

- The PATIENT is an intended OPERATOR.
- The device must only be serviced and maintained when switched off.
- The device may be maintained by the user itself as described in the according section of this manual.
- Stop using the equipment immediately after contact with water.

REGARDING EMC STATEMENT

1. A statement to the effect that this equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS; and
2. a statement to the effect that wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept distanced away from the equipment. The distance is calculated by the MANUFACTURER from the 800 MHz to 2,5 GHz column of Table 5 or Table 6 of EN 60601-1-2:2015, as appropriate.

MANUFACTURER'S WARRANTY

The producer warrants Wellion WAVE blood pressure meter against defects in materials and workmanship for a period of two years from the date of purchase. The warranty is lost if the instrument is misused, poorly maintained or is tampered with.

Liability under this warranty is limited to the repair of defective parts or - at the discretion of the producer - to the replacement of the instrument. The right to rescind the purchase agreement exists only if the replacement is also found to be defective. Claims other than these will not be dealt with.

The warranty is not valid if the damage results from misuse, maltreatment, tampering, human error and the use of extreme force. This warranty is valid only if the date, the stamp and the signature of the dealer are recorded on the warranty card on the date of purchase. The warranty period is not to be extended by any claim made under this warranty.

Appendix 1 Guidance and Manufacturer Declaration Tables

Guidance and manufacturer's declaration- electromagnetic emission		
The Model PG-800A11 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A11 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Model PG-800A11 Series Electronic Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complied	The Model PG-800A11 Series Electronic Blood Pressure Monitor is used in home and it's powered by DC 3V.

Guidance and manufacturer's declaration – electromagnetic immunity			
The Model PG-800A11 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A11 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m,50/60Hz	30 A/m,50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a. c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration-electromagnetic immunity

The Model PG-800A11 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A11 Series Electronic Blood Pressure Monitor should assure that it is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 6100-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bandsa	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the Models PG-800A11 Series Electronic Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{V_t} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 80MHz to 800MHz 800MHz to 2.7GHz
Radiated RF IEC 6100-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b Interference may occur in the vicinity of equipment marked with the following symbol: ((⋯)) 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model PG-800A11 Series Electronic Blood Pressure Monitor is used exceeds the applicable RF compliance level above, the Model PG-800A11 Series Electronic Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model PG-800A11 Series Electronic Blood Pressure Monitor.

d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Model PG-800A11 Series Electronic Blood Pressure Monitor

The Model PG-800A11 Series Electronic Blood Pressure Monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model PG-800A11 Series Electronic Blood Pressure Monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model PG-800A11 Series Electronic Blood Pressure Monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{12}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.2	1.2	1.2	2.3
10	3.8	3.8	3.8	7.3
100	12	12	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

CALIBRATION METHOD

1. Press and hold the „ON/OFF and M“ button at the same time, load the battery, enter the static air pressure calibration mode after the LCD screen is fully displayed, and then release the button.
2. Press ON/OFF to code the internal air valve.
3. Connect the external standard barometric interface and the digital barometer interface to the cuff interface.
4. External input 50mmHg and 200 mmHg standard static air pressure, and observe the air pressure value displayed at the position of the LCD systolic pressure (SYS) and the value of the digital pressure gauge should be in the range of ± 3 mmHg.

Caution:

1. ME devices can be used in exposed environments, including electromagnetic interference environment to ensure basic safety and basic performance unchanged.
2. In the event of any serious event related to this product, such as serious adverse event, significant alteration of the product resulting in change of intended use, etc., it will be reported to the manufacturer and the competent authorities of the user and/or the member states where the patient is located.

Notes:

Essential performance: Limits of the error of the manometer, ± 3 mmHg. Reproducibility of the blood pressure determination ± 3 mmHg.

Clinical benefits: Accurate measurement of SBP and DBP, clinical performance meets the requirements of ISO 81060-2:2018.

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© MED TRUST, Wellion and WAVE are registered trademarks.